

The world's first [8D Anti-aging HIFU]
The overall effect treatment for face
lifting and skin tightening and
contouring.
Zero downtime

Single row + Multi row + Multifrequency + 360°line carving



Cartridge



FL Handle



SD Handle



Single row handle working head
(optional for S2.0)



Multi row handle working head
(optional for WD5.5, WD9)



R1.5mm




R3.0mm



R4.5mm

(optional Configuration)



Dr Genesis 8D HIFU mechanism actions



Advantage:

Global voltage mature and stable technology! Depth precision, energy intelligent adjustment; intelligent self-checking system to ensure the best output state. The state presented from the outside to the inside can be said to be a very good one beauty instrument.

For loose skin and want to regain firmness!

Facial Firming - Powerful tightening, lifting eyebrows, nasolabial folds

Body Firming - Lifting to tighten arms, stomach, thighs

Improve skin elasticity - cell regeneration, elasticity enhancement


Facial single row precision + 360 ring output, according to the facial special micro and no dead angle lifting, energy penetrates deep into the skin, stimulates collagen regeneration, makes cells self-repair, and finally achieves the effect of wrinkle rejuvenation, especially for lifting facial skin. It has obvious effect on diluting nasolabial folds and fine lines at the bottom of the muscle. In addition, it stimulates the skin to restore cell activity, enhances blood circulation, and eliminates dark circles. Burn excess fat. Multiple rows for each part of the body, efficient and precise.

Depth precision CFU 360 ring output

Biotechnology for extreme lifting

What is CENTERLESS FOCUSED ULTRASOUNDS (CFU)

- The Centerless Focused Ultrasound creates a Thermal Coagulation Point (TCP) in the dermis layer for a non-invasive and safe lifting effect
- The Centerless Focused Ultrasound emitted the heats precisely into the point of focus, where it reaches a temperature to 65-70°C to stimulate new collagen and elastin fibers without damaging the surrounding skin.
- The Centerless Focused Ultrasound safely and effectively creates a perfect oval TCP. According to experimental results, it has shown to reduce side effects as well as improving the performance of the heat lesion by more than 50%.




Normal HIFU

8D HIFU

Part 1 – Skin Layer

The high intensity focused ultrasound energy penetrates through the temperature rise above 65°C to the treated area, at the desired depth it generating a thermal coagulation effect.

It is non-invasive treatment that helps stimulates collagen and elastin fibers without damaging the surrounding tissue, in addition it also activate the contraction of muscle fibers providing a perfect V-shape and giving the skin plump, firm and lift



8 Depth of Dimension

